

Beijing Huacheng Taike Technology Co., Ltd

Web:www.cosmedplus.com

teklaser.en.alibaba.com

E-mail:info@cosmedplus.com

Tel:0086(010)-61592593

Phone:0086 18701668578

Fax: 0086(010)-61592593

Cosmed Plus Cmp GmbH

Kreissparkasse Heidenheim Schnaitheimer

Str. 12 89520 Heidenheim

USt-IdNr. DE 145 564 690

Tel:07321/344-0



C4
COSMEDPLUS



PRODUCT MANUAL



BEIJING HUACHENG TAIKE TECHNOLOGY CO.,LTD

”” BUSINESS PROFILE

Established in 2005, Beijing Huacheng Taike Technology Co., Ltd. is a professional manufacturer of beauty and medical equipment. We have been specializing in the research and development, manufacturing, marketing and service of beauty product lines for more than ten years and gained rich experience. In 2017, we set up a Research and Development Center in Germany, aiming to take full advantage of the world's most advanced technology and R&D capabilities to bring the company's innovation to a higher level. We are focusing on more efficient, more convenient service for European customers, better to meet customer needs.

So far, we have won thousands of customers in the UK, Poland, Czech Republic, Germany, Italy, Turkey, Iran and other countries and regions. Meanwhile, we provide OEM services to some of these exclusive agents.

Huacheng Taike is committed to the highest level of customer service, competitive prices, fast delivery and comprehensive, cutting-edge products. To win your satisfaction is our ultimate goal.



Contents

- 1-2 Portable Hydra facial
- 3-4 Vertical Hydra facial
- 5-6 1200w+600w diode laser hair removal machine
- 7-8 600w diode laser hair removal machine
- 9-10 1200w diode laser hair removal machine
- 11-12 2400w diode laser hair removal machine
- 13-14 Cryotherapy
- 15-16 Vertical Q-switch
- 17-18 Portable Q-switch
- 19-24 4D HIFU High Intensity Focused Ultrasound
- 25-26 Home use hair removal machine
- 27-28 Co2 laser
- 29-30 Pico second
- 31-32 IPL
- 33-34 CMP-CL series air cooling machine

Function

Deep skin clean.
 Detox.
 Facial skin remodeling and improve skin elasticity.
 Anti -Aging.
 Elimination of wrinkle.
 Skin rejuvenation.
 White / Black heads removal, deep cleaning
 Moisturizing
 Whitening , Lightening and Skin Rejuvenation
 Skin Reproduction
 Keratin Removal

Handle 1:Radio frequency skin lifting
 (RF Skin tightening, Stimulate collagen)
 Handle 2:Cold Hammer-skin improve
 (Skin Rejuvenation, shrink skin pores)



Handle 3:Dead skin Removal-Skin Repair
 (Remove dead skin, improve Skin)
 Handle 4:H2 O2 skin Deep Cleaning
 (Remove black head, remove skin grease dirty, improve oily skin, deep skin cleaning,increase skin elasticity, improve complexion)

Handle 5:H2 O2 Skin Moisture
 (Replenish skin oxygen and water skin rejuvenation, fast inject nutrient to skin increase cell activity)
 Handle 6:Ultrasonic-Nutrient Inject
 (Skin deep permeate,Skin repair)
 Skin Scrubber handle
 Use high frequency vibration wave, softens the skin cutin cells to clean the dirty deeply
 Hydro Facial handle
 Clean the cuticle and hair follicle layer to make the skin shiny
 Multi-RF handle
 RF made skin absorb lots of heat energy, Produce new collagen, smoothing wrinkles and tightening the skin
 Cold hammer
 Use for clam down and repair damaged skin ,eliminate redness and sensitive effect
 Ultrasonic handle
 Ultrasound vibrate 1 million to 3 million times per Minute to opening pores,the essence of rapaid infiltration to the base of skin, skin absorption rate reaching more than 90%

Specification

Tornado tip	Large, Small
Vacuum range	Max 680 mm
Radio frequency	1Mhz, Multi polar(6Pol)
Lon lifting	500Hz (Digital Lifting)
Ultrasound	1Mhz / 2W/cm2
User interface	7inch Color Touch LCD
Input voltage	100~120V, 50/60HZ / 220~240V, 50/60HZ
Noise level	45Db
Weight	12KG



7 in1 Vertical Hydra facial Pro

Applications:

- Deep skin clean
- Detox
- Facial skin remodeling and improve skin elasticity
- Anti –Aging
- Elimination of wrinkle
- Skin rejuvenation

Features:

Handle 1 Water Spray Gun

The instrument generates high -oxygen -water in combine solution, and the spray gun lead into the dermis layer.



Handle 2 Two-pole RF

Promote collagen regeneration fill wrinkles, restore elasticity.

Handle 3 Aqua peeling head

Pure water is purified to produce H2 fine molecules on the surface of the skin, so what water molecules quickly penetrate the skin

Handle 4. Hot and cold desensitization handle

Select the cooling or heating function to adjust the temperature. Work the front part of the handle to the active part.

Handle 5 Multipole RF

Using heat to stimulate the dermis to secrete more collagen to fill the shrinking collagen vacancies, hold up the skin scaffold, restore skin elasticity, make the skin rich, fresh and vibrant.

Handle 6 Carbon Oxygen

The skin is in a clean and very moist state, using a membrane on the carbon-oxygen kinetic energy massage head combined with bright white muscle gel and vitality regenerating gel to the skin, the product and the membrane combine to produce rich small bubbles, forming a Boer effect. This is the oxygen kinetic energy technology. The epidermis produces a large number of small bubbles, that is, carbon dioxide, which causes the skin to be in a weak acid environment. At this time, the pH value in the cells decreases, the affinity of hemoglobin oxygen in the red blood cells decreases, and hemoglobin releases more oxygen, promotes facial blood circulation, and accelerates metabolism.

Handle 7 Ultrasonic water light introduction---nutrition reserve

Through the introduction of needle-free water, different types of nutrients are introduced into the deep skin to awaken the cells and repair the withered and broken cells, thereby strengthening the skin lock water. Reshapes firm, smooth, elastic skin in a short period of time.

Vital Firming Treatment → Skin Regeneration Treatment

Brightening skin care → even whitening skin tone





New generation diode laser device

Introduction

The 755nm+808nm+1064nm diode laser represents a new frontier in laser technology for permanent hair removal. This device creates a laser beam with a wavelength of 755nm+808nm+1064nm that is able to pass through the skin and gradually cut off the blood supply to the hair bulb as it is absorbed by the pigment in the hair and raises the temperature. So that hair will disappear without injury surrounding tissue.



Application

Permanent hair removal

Technical Advantages

1. America Coherent laser, the best laser manufacture in the world which guarantee 40,000,000 shots
2. Japan TEC cooling system ,makes the cooling faster and more efficiency
3. Water filter secured to ensure long life of the whole system.
4. Modularized design for convenient maintenance.
5. 15.6 inch colorful android touch screen meets the Age of Intelligence
6. Sturdy injection molded shell.
7. Appearance design team from Italy makes machine appearance cater to the world's aesthetic standards
8. Important parts are manufactured overseas and certified by medical CE.



Specification

Wavelength	755nm+808nm+1064nm
Laser Output	600W/1200W
Frequency	1-10Hz
Spot Size	2cm2/4cm2
Pulse Duration	1-300ms/1-400ms
Energy	1-120J /1-180J
Cooling System	Japan TEC cooling system
Sapphire contact cooling	-5-0℃
Operate Interface	15.6 inch color touch android screen



New generation diode laser device

Introduction

The 755nm+808nm+1064nm diode laser represents a new frontier in laser technology for permanent hair removal. This device creates a laser beam with a wavelength of 755nm+808nm+1064nm that is able to pass through the skin and gradually cut off the blood supply to the hair bulb as it is absorbed by the pigment in the hair and raises the temperature. So that hair will disappear without injury surrounding tissue.


Application

Permanent hair removal

Technical Advantages

1. America Coherent laser, the best laser manufacture in the world which guarantee 40,000,000 shots

2. Japan TEC cooling system, makes the cooling faster and more efficiency
3. Water filter secured to ensure long life of the whole system.
4. Modularized design for convenient maintenance.
5. 15.6 inch colorful android touch screen meets the Age of Intelligence
6. Sturdy injection molded shell.
7. Appearance design team from Italy makes machine appearance cater to the world's aesthetic standards
8. Important parts are manufactured overseas and certified by medical CE.



Specification

Wavelength	755nm+808nm+1064nm
Laser Output	600W
Frequency	1-10Hz
Spot Size	2cm ²
Pulse Duration	1-300ms
Energy	1-120J
Cooling System	Japan TEC cooling system
Sapphire contact cooling	-5-0°C
Operate Interface	15.6 inch color touch android screen



New generation diode laser device

Introduction

The 755nm+808nm+1064nm diode laser represents a new frontier in laser technology for permanent hair removal. This device creates a laser beam with a wavelength of 755nm+808nm+1064nm that is able to pass through the skin and gradually cut off the blood supply to the hair bulb as it is absorbed by the pigment in the hair and raises the temperature. So that hair will disappear without injury surrounding tissue.

Application

Permanent hair removal

Technical Advantages

1. America Coherent laser, the best laser manufacture in the world which guarantee 40,000,000 shots
2. Japan TEC cooling system,makes the cooling faster and more efficiency

3. Water filter secured to ensure long life of the whole system.
4. Modularized design for convenient maintenance.
- 5.15.6 inch colorful android touch screen meets the Age of Intelligence
6. Sturdy injection molded shell.
7. Appearance design team from Italy makes machine appearance cater to the world's aesthetic standards
8. Important parts are manufactured overseas and certified by medical CE.



Specification

Wavelength	755nm+808nm+1064nm
Laser Output	1200W
Frequency	1-10Hz
Spot Size	4cm ²
Pulse Duration	1-400ms
Energy	1-180J
Cooling System	Japan TEC cooling system
Sapphire contact cooling	-5-0°C
Operate Interface	15.6 inch color touch android screen



New generation diode laser device

Introduction

The 808nm+ 1064nm diode laser represents a new frontier in laser technology for permanent hair removal. This device creates a laser beam with a wavelength of 808nm+ 1064nm that is able to pass through the skin and gradually cut off the blood supply to the hair bulb as it is absorbed by the pigment in the hair and raises the temperature. So that hair will disappear without injury surrounding tissue.

Application

Permanent hair removal

Technical Advantages

1. America Coherent laser, the best laser manufacture in the world which guarantee 40,000,000 shots
2. Japan TEC cooling system which makes the cooling faster and more efficiency

3. Water filter secured to ensure long life of the whole system.
4. Modularized design for convenient maintenance.
5. 15.6 inch colorful android touch screen meets the Age of Intelligence
6. Sturdy injection molded shell.
7. Appearance design team from Italy makes machine appearance cater to the world's aesthetic standards
8. Important aparts are manufactured overseas and certified by medical CE.



Specification

Wavelength	808nm/808+1064nm/755+808+1064nm
Laser Output	2400W
Frequency	1-3Hz
Spot Size	22*35
Pulse Duration	1-300ms/1-400ms
Energy	5-60J
Cooling System	Japan TEC cooling system
Sapphire contact cooling	-5-0°C
Operate Interface	15.6 inch color touch android screen
vacuum assisted	0.03-0.09Mpa



Introduction

Cryolipolysis is a method of placing a cryolipolysis instrument on the surface of human skin to cool subcutaneous tissues, premature aging of fat cells, and death in succession. It is excreted out of the body through metabolism to achieve a slimming effect.

The process of freeze fat dissolving gradually absorbs the heat of subcutaneous fat. The temperature of fat cells is reduced to zero degrees Fahrenheit, which means that they are frozen. Low temperature does not affect skin or muscle while killing fat cells. Dead fat cells are then excreted through the liver.

For those who are overwhelmed with "stubborn" fat, cryolipolysis is undoubtedly a godsend. Whether for fat-intensive parts or parts with a small fat area, such as love muscles (loose fat on the sides of the waist above the hips), belly fat, and back fat, this weight loss technique can produce stunning results.

The Applications of Cryotherapy

Cellulite Loss, Freeze fat , Cryotherapy

The Advantages of Cryotherapy

- 1.The current non-surgical medical cosmetic technologies, each treatment, patient can lose 5- 8kg.
- 2.Four different sized fat freeze handpieces,can be used for big areas and small areas.
- 3.More advanced than liposuction, fat melting technology.
- 4.Big water tank, cooling is very good.
- 5.All handpieces are water-proof design, safe and stable.
- 6.The freezing suction cups are made of silicone, soft and comfortable when touching the skin.

The Specifications of Cryotherapy

Product name	Loss weight machine equipment
input voltage	AC110V/220V 50-60Hz
power consumption	<500w
screen	10-ch wide color touch screen
cooling device output temperature	-11-45°C
cooling device output pressure	0-100KPa
cooling liquid	pure water
default time	0-120min
fuse size	F2AL250V
Ambient temperature	5°C-+40°C
relative humidity	≤80%
carton size	60x57x138cm
atmospheric pressure	80KPa-106KPa
G/w	95KG
warranty time	one year
delivery time	7-14 days
terms of payment	T/T with 30% deposit in advance, the balance should pay before shipment
Patent NO.	ZL201130025860.8
MOQ	1 set
power	Less than 450VA



Vertical q-switch

Introduction

Utilizing the explosive effect of the ND: YAG LASER, the laser permeate the By epidermis into the dermis and takes effect on pigment mass and is absorbed pigment. Since the laser pulses is extremely short in nanosecond and comes with super high energy, the pigment mass would swell quickly and break into small pieces, which will be eliminated through metabolism. Thus the tattoo or other pigmentation will be removed without injury to normal tissue. The treatment is safe and convenient without downtime or side effects.



Applications

1. Tattoo Removal (including eyebrow tattoo, eyeliner, lip line)
2. Pigmentation Removal
3. Skin Rejuvenation

Technical Advantages

1. USA lamp which can guarantee 2,000,000 shots
2. Infrared indication makes more professional
3. Multi languages supported, meets requirements of global market
4. Separate water system and electronic circuit, safe for transportation
5. Auto-alert for low water level, water temperature and water flow.
6. Plug and play handle.

Specification

Power	500w/1000W
Wavelength	532nm/755nm/1064nm/1320nm
Energy	0-2000mj
Frequency	1-10HZ
Pulse duration	<6ns
Yag bar	Φ7
Operate interface	10.4 inch color touch screen
Cooling system	Air cooling & water cooling
Power supply	AC85-130V or AV180-260V, 50/60HZ



Portable q-switch

Introduction

Utilizing the explosive effect of the ND: YAG LASER, the laser permeate the By epidermis into the dermis and takes effect on pigment mass and is absorbed pigment. Since the laser pulses is extremely short in nanosecond and comes with super high energy, the pigment mass would swell quickly and break into small pieces, which will be eliminated through metabolism. Thus the tattoo or other pigmentation will be removed without injury to normal tissue. The treatment is safe and convenient without downtime or side effects.



Applications

1. Tattoo Removal(including eyebrow tattoo, eyeliner, lip line)
2. Pigmentation Removal
3. Skin Rejuvenation

Technical Advantages

1. USA lamp which can guarantee 2,000,000 shots
2. Infrared indication makes more professional
3. Multi languages supported, meets requirements of global market
4. Separate water system and electronic circuit, safe for transportation
5. Auto-alert for low water level, water temperature and water flow.
6. Plug and play handle.

Specification

Power	500w
Wavelength	532nm/755nm/1064nm/1320nm
Energy	0-2000mj
Frequency	1-10HZ
Pulse duration	<6ns
Yag bar	Φ6
Operate interface	8 inch color touch screen
Cooling system	Air cooling & water cooling
Power supply	AC85-130V or AV180-260V, 50/60HZ



HIFU (High Intensity Focused Ultrasound) is the most effective method currently available for the renewal of collagen, one of the latest advances in the beauty industry.

On HIFU technology, the ultrasonic waves emitted by the treatment head intersect at one point, so their combined effect is intensified.

During treatment with HIFU, this microfocused ultrasound reaches the skin's various layers and induces a coagulation process at 65-75°C temperature, triggering collagen regeneration in tissues.



What makes 3rd Generation HIFU more efficient than traditional equipment?

Developers of 3rd Generation HIFU have managed to achieve more accurate results than ever before, with the accuracy of the software calculated by the operator. As shown in the figure below, it is due to this precision that collagen can evenly build up on such a large surface. With the 3rd Generation HIFU, the operator can handle the treatment much more readily than the conventional 1-line operator operating machines.

For conventional techniques, not only does the desired result fall, but the undesirable effect can be caused by the superficial treatment resulting from the irregular layout. With the 3-generation HIFU professional treating head and precision software, the 10-line arrangement allows collagen fibers to form evenly in the deep tissue of the skin tissue, which naturally also exfits the skin evenly to provide maximum results with HIFU treatment. Additionally, the 10-line treatment head can handle a larger area for one treatment, thus increasing the effectiveness of the treatment.


Face Lifting



This treatment targets the upper layers of the skin and dermis for collagen regeneration, via advanced ultrasound energy at the depth of 1.5mm. This is perfect for use around eye and lip areas. This treatment offer results similar to fractionated needling device services, but without the use of needles for skin penetration, avoiding any downtime.

Body Slimming





Benefits

- Refines lines and wrinkles
- Tones skin and reduces the size of pores.
- Reduces the appearance of Acne scarring
- Reduces the appearance of stretch marks
- Best results are achieved with a series of 4-6 sessions at 6 week intervals.

HIFU Facelift treatment

The closest treatment to a surgical facelift. High frequency Ultrasound beams pass through the skin, to deliver precise and fractional ultrasound energy, at exactly the right depths and temperature, to stimulate both muscle firming/lifting and new collagen production, all with no damage or disruption to the skin's surface. Consider this treatment a natural facelift. It works by the connective tissue repairing itself after be treated with the ultrasound beam.

Body Slimming




Looking for HIFU fat removal?

HIFU offers a real natural alternative to liposuction or surgery. Only one treatment is required to achieve the results with no downtime. HIFU body sculpting is a non-invasive treatment for removing stubborn body fat, skin toning and tightening.

Treatment areas

- Under arms (Tuck shop arms)
- Buttocks
- Inner/outer thighs
- Thighs
- Stomach
- Perfect for localised fatty deposits



HIFU (High Intensity Focused Ultrasound) penetrates through the surface of the skin to deliver precise and fractional focused ultrasound energy, into the fat layer of the skin (6.0 & 8.0mm – underarm / back areas and 13.0 & 16.0 mm for the abdomen, buttocks, thigh areas), and at the right temperature (60~70°C), without damaging or disrupting the epidermal tissue (surface) of the skin. HIFU (High Intensity Focused Ultrasound) penetrates through the surface of the skin to deliver precise and fractional focused ultrasound energy, into the fat layer of the skin (6.0 & 8.0mm – underarm / back areas and 13.0 & 16.0 mm for the abdomen, buttocks, thigh areas), and at the right temperature (60~70°C), without damaging or disrupting the epidermal tissue (surface) of the skin.

Functions

- Remove wrinkles on around forehead, eyes, mouth, etc.
- Lifting and tightening both cheeks skin.
- Improving skin elasticity and shaping contour.
- Tightening the skin tissue on forehead, lifting the eyebrows lines.
- Improving skin complexion, making the skin delicate and bright
- Removing neck wrinkles, protecting neck aging.



HIFU Cartridges arrange	HIFU Cartridges types	Treatment areas
Face lift cartridges	DL-1.5mm	Eyes around
	DL-3.0mm	Cheek/Neck
	DL-4.5mm	Chin
Body slimming cartridges	DL-6.0mm	Chest/Arms
	DL-8.0mm	Arms/legs
	DL-10.0mm	Abdomen/Waist/Buttock/Thigh
	DL-13.0mm	Abdomen/Waist/Buttock/Thigh
	DL-16.0mm	Abdomen/Waist/Buttock/Thigh



Specifications

Energy output	0.2-8.0J(0.1J/step)
Spacing	1.0-10mm(0.5mm/step)
Length	5.0-25mm(1.0mm/step)
Cartridge optional	1.5mm,3.0mm,4.5mm, 8mm,11mm,13mm,16mm (standard 3 cartridge)
Cartridge lift time	No limit shots

808nm home use laser hair removal device

Performce parameter

- 1.wavelength:808nm
- 2.laser cooling system:forced air cooling
- 3.spot size:12mm*15mm,sapphire crystal,single-side coating
- 4.laser power:0-50w(0-12J)/0-100w(0-24J)
- 5.display interface:1.8 inch colorful screen
- 6.adapter inter voltage:AC180~270V
- 7.adapter power:240w
- 8.use environment temperature:-10℃ ~+40℃
- 9.save environment temperature:-40℃ ~+60℃
- 10.light outlet ice compress function
(cooling power 10w,temperature 15℃ ~20℃)
- 11.ice compress texture:aluminium alloy
- 12.shell surface:mirror spray paint
- 13.packaging:flight case
- 14.box size:420X300X170(mm)
- 15.product weigh(including packaging and parts):less than 10kg



Graphical representation

Functional description

- 1.Ice compress head with alarm for high temperature and temperature detection.
- 2.Laser temperature detection with power and over temperature protect.
- 3.Body contact induction and light control switch make it safer.Only when the device touch your skin and its switch is pressed,it will work.
- 4.Long time and ongoing light.




Gear control

100w 2bar

Level	Power	Pulse length
1st	10w	80ms
2nd	20w	160ms
3rd	30w	240ms
4th	40w	320ms
5th	70w	400ms

Instructions


- 1.push the device switch over 2s to turn on the device until the **Welcome** is appear.Then 2s later **808** is appear and enter the site selection interface
 - 2.select the site
 - 3.you can turn up or down the power with button
- Facial interface
- Body interface
- 4.push the device switch over 2s to turn off





CO2 Fractional Laser

The new developed RF-Excited CO2 fractional laser fires a laser beam through RF CO2 tube which is split into numbers of microscopic beam, producing tinier dot than general CO2 laser (glass tube). The treatment tip can vaporize the outmost layers of the skin across a large surface by creating thousands of tiny microscopic laser wounds, spaced in an even pattern across the skin, but leaving areas of healthy, untreated skin between them, the lower collagen layer of the dermis is stimulated to renew and repair. Therefore, the heat of laser only passes deeply through the factional damaged area, the surface of skin now contains only microscopic surface wounds, in stead of a large, red, oozing burn. During the skin self-resurfacing, a great amount of collagen is produced for skin rejuvenation, after some recovery, the newly-generated skin is remarkably smoother, healthier.



Clinical Application

1. Smooth scars such as surgical scars, burnt scars, acne scars, etc
2. Skin renewing and resurfacing, sun damage recovery
3. Wrinkle removal and skin tightening
4. Pigmentation removal such as intractable chloasmas, age spots, blemishes, speckle etc
5. Skin excrescence such as wart,etc
6. Acne treatment



Technical Advantages

1. The self-test system of the machine will immediately start up after turn it on; the self-test
2. System will be continuous monitoring of the circuit during equipment operation. The laser transmit to the patient's skin via fractional head, the laser array mirror only shot from the front face handle.
3. The machine has designed a separate safety circuit, when there is leakage of the circuit ,it can be cut power off automatically.
4. Emergency appears that can be quickly cut off the power by the red emergency switch.
5. Key switch to prevent any person free to open the device.


Co2 laser specification

Laser type	RF Tube(USA Coherent)
Output power	30w
Wavelength	10600nm, far-infrared laser
Pulse energy	1mj-400mj
Scan graphics(7graphics)	Square,rectangle,round,triangle,oval,6-diamond shape,line
Scan graphic size	0.1*0.1mm-20*20mm
Scan mode(3 modes)	Random,sequence,middle to sides
Scan interval	0.1-6s
Spot distance	0.1-2.0mm
Pulse width	0.1-10ms
Aiming beam	635nm/10mw
Power supply	220v/110v
Dimension	58*55*126cm
Net weight	45kg



application

The Q-Switched (QS) laser is the most effective way to remove natural or artificial (tattoo) pigmentation, while minimizing the risk of damage to surrounding tissue. Uses an active Q-switched laser to deliver photo acoustic shockwaves to the target area through high laser intensities in short nanosecond pulses. The mechanical Q-switched effect works by vibrating and breaking up the pigment in the lesion or the ink particles in the tattoo. As the area heals, the body's immune system flushes away the shattered pigment, revealing lighter, clear skin with minimal risk of scarring or hyper pigmentation.



Technical advantages

- 1.Pigment clear rate is higher
- 2.Pigment removal speed is faster
- 3.Refused to stimulate melanin active again risk
- 4.Remove melanin colleagues and staart repair mechanism
- 5.Don't damage the normal tissue



Specofocatopn

SINGLE YAG BAR Φ6/Φ7

Power	750W
Wavelength	532nm;1064nm
Energy	0-2000mj
Frequency	1-20HZ
Pulse duration	500PS
Yag bar	Φ6/Φ7
Operation interface	10.4 inch color touch screen
Spot size	2mm-10mm
Power supply	≥750W
Dimension	117*54*121(cm)
Net weight	76KG

DOUBLE YAG BAR Φ6&Φ7

Power	1200W
Wavelength	532nm;1064nm
Energy	0-2000mj
Frequency	1-20HZ
Pulse duration	500PS
Yag bar	Φ6&Φ7
Operation interface	10.4 inch color touch screen
Spot size	2mm-10mm
Power supply	≥1200W
Dimension	117*54*121(cm)
Net weight	80KG



Before After Before After Before After



Vertical OPT SHR E-light Device

Introduction

IPL is a kind of light with high energy and broad wavelength, with the range of 420nm to 1200nm. Intense pulsed light can act directly on the skin, make skin lesions on the solar energy selective absorption to the treatment of a variety of tissues, and do not damage normal skin tissue.

Introducing Super Hair Removal Technology (SHR) Lasers - a revolutionary new method for permanent hair removal that is virtually pain & side effect free. In comparison to other slightly dated laser and IPL methods, SHR provides faster, safer and pain free hair reduction treatments for clients.



Technical Advantages

1. UK lamp which guarantee 300,000 shots.
2. Advanced cooling system, adjustable contact cooling for comfortable treatment.
3. Japan Nippon capacitors
4. German medical power makes the energy stable.
5. Multi-languages supported, meets requirements of global market.
6. Sturdy injection molded shell.
7. Friendly system with standard mode and professional mode.



Application

1. Hair removal
2. Acne removal
3. Wrinkle removal
4. Freckle removal
5. Vascular removal
6. Pigment removal

Specification

Wavelength	430/530/560/610/640~1200nm(standard) 690~950nm for SHR mode 480/585/690/755~1200nm(Optional)
Spot size	15*50mm
IPL light Energy	10~60J/cm2
Pulse Duration	2~20ms
Pulse Sequence	1~5 pulses
Delay time	5~60ms
Intergrated skin cooling	-4~0℃
Operation interface	8" ture color touch screen
Drive Power	2,400W with SHR mode
Working period Continuously	8~10 hours
Power supply	AC85~130V or AC180~260V, 50/60HZ
Cooling system	Air&water cooling circulation with copper radiator treatment head with semiconductor cooling
Weight kgs	49kg



CMP-CL Series air cooling machine - low noise, strong comfort

Introduction: Compressing air to compress the air and then injection cold air through the refrigerating device. Continuously provide cold protection to the epidermis, protect the epidermis from thermal damage caused by photothermal action, to protect the epidermis, reduce pain and not damage the epidermis. The purpose of making high energy density on the target tissue.

Features: Almost all phototherapy and cosmetic surgery require good protection of the epidermis. CMP-CL Series air cooling machine by providing -10°C .The following cold air can provide cold protection to the epidermis before, during and after surgery, so that the epidermis is protected from thermal damage caused by photothermal action. The thermal damage that has been generated can also be mitigated by the cold protection provided.

Advantages

1. After using Co2 or Q-Switch laser, use air-cooled machine to reduce redness and make the process more effective and smooth.
2. For thermal damage caused by OPT rapid hair removal treatment, an air cooler can be used to reduce skin temperature.

Application: This product is suitable for all kinds of skin treatment instruments, beauty instruments, weight loss instruments, large photon instruments, tattoo removal machines, eyebrow removal machines, etc. It is widely used as an auxiliary instrument in the medical and cosmetic industries.



Specification

Model	CMP-CL10	CMP-CL15	CMP-CL25 CMP-CL25A
Power supply	220/110V	220/110V	220/110V
Air volume	0-2m3/min No shifting	0-2m3/mi No shifting	0-2m3/min 9-speed shift
Power	800W	1000W	1500W
Noise	<60db	<60db	<60db
Refrigerant	404A Environmentally friendly	404A Environmentally friendly	404 Environmentally friendly
Lowest temperature	-15°C	-25°C	-35°C
Weight	30kg	50kg	75kg
Dimension	350*300*450mm	310*550*700mm	430*550*780mm
Package dimensions	450*400*550mm	400*650*800mm	560*770*860mm