



#### Vacuum 3D HIFU

- 1.more comfortable feeling with vacuum suction
- 2.treatment lines of one press;1 line-11 lines;
- 3.cartridge; 3.0mm,4.5mm,6.0mm,8.0mm, 10.0mm,13.0mm,16.0mm treatment depth;
- 4.the energy shots be more well-distributed,so treatment results will be better than others;
- 5.total shots;10000lines and 500lines(spare);
- 6.anti-aging,wrinkles removal,lifting and tightening skin



#### 2D HIFU

- 1.treatment width of one step:1mm-10mm;
- 2.treatment lines of one press:1line-11lines;
- 3.cartridge: 3.0mm , 4.5mm ,6.0mm,8.0mm, 10.0mm,13.0mm,16.0mm treatment depth ;
- 4.the energy shots be more well-distributed,so treatment results will be better than others;
- 5.total shots:10000lines and 500lines(spare);
- 6.anti-aging,wrinkles removal,lifting and tightening skin.




#### 2D HIFU (ABS)

- 1.treatment width of one step:1mm-10mm;
- 2.treatment lines of one press:1line-11lines;
- 3.cartridge: 3.0mm , 4.5mm ,6.0mm,8.0mm, 10.0mm,13.0mm,16.0mm treatment depth ;
- 4.the energy shots be more well-distributed,so treatment results will be better than others;
- 5.total shots:10000lines and 500lines(spare);
- 6.anti-aging,wrinkles removal,lifting and tightening skin.



#### HIFU for treating vagina

- 1.360° rotation automatically;
- 2.cartridge: 3.0mm , 4.5mm treatment depth ;
- 3.strengthen concentrating vagina's forces,
- 4.promoting collagen reborn;
- 5.tightening vagina;
- 6.Increasing vaginal elasticity;
- 7.moisturizing vagina;
- 8.improving incontinentia urinae.



#### Vagina HIFU

- 1.360° rotation automatically;
- 2.cartridge: 3.0mm , 4.5mm treatment depth ;
- 3.strengthen concentrating vagina's forces,
- 4.promoting collagen reborn;
- 5.tightening vagina;
- 6.Increasing vaginal elasticity;
- 7.moisturizing vagina;
- 8.improving incontinentia urinae.

## CATALOGUE

Vacuum 3D HIFU

2D HIFU

2D HIFU (ABS)

HIFU for treating vagina

Vagina HIFU

Thermolift Focus rf

Liposonix HIFU Vacuum rf machine

VEVIVE

Thermolift with vacuum

Liposonix HIFU with RF and vacuum

M5

Water and oxygen beauty equipment

Qo

BIO Micro-current

Remage (thermage)

Fly Action

Dream Action Plus(portable)

Qo

Water and oxygen beauty equipment

# MENO®

MENOBauty Technology Co.,LTD.

www.menobeauty.net  
www.menobeauty.cn

Add:6F,XinBao Yi trade and industry building,NO32 the 3rd industrial,HouTing Village,AhJing Town,BaoAn District,ShenZhen 518000,China

Tel : ( +86 ) 0755-27838931/27883993

Fax : ( +86 ) 0755-27838931/27883993Ext601

# MENO®

MENOBauty Technology Co.,LTD.



PRODUCTS DESIGNED AND PRODUCED CAREFULLY,LEADER IN THE BEAUTY LINES

# HIGH-END BEAUTY MACHINES

## HIGH-END BEAUTY MACHINES

## MANY NEW MACHINES ARE COMING



#### Thermolift Focus rf

- 1.40.68MHZ rf frequency unipolar rf+bipolar rf (unipolar rf+vacuum rf)Unipolar rf;
- 2.thermotherapy from SMAS to epidermis, power 300w Bipolar rf thermotherapy from epidermis to SMAS ,power 300w Skin lifting, tightening and body shaping,weight loss.



#### Liposonix HIFU Vacuum rf machine

- 1.40.68MHZ rf frequency unipolar rf+bipolar rf (unipolar rf+vacuum rf)Unipolar rf;
- 2.thermotherapy from SMAS to epidermis, power 300w Bipolar rf thermotherapy from epidermis to SMAS ,power 300w Skin lifting, tightening and body shaping,weight loss.



#### VEVIVE

- 1.Thermolift Focus RF technology with 40.68MHZ frequency;
- 2.Non-surgical,painless medical ;
- 3.treatment for treating vaginal laxity;
- 4.Works by cooling and gently heating vaginal;
- 5.tissue and increases feeling of vaginal tightness and sexual satisfaction;
- 6.Treatment is completed in 30 min,no downtime.




#### Thermolift with vacuum

- 40.68MHZ unipolar radio frequency and vacuum; lift face skin,improve skin color, tighten and shape body;5pcs tips for face or body.



#### Liposonix HIFU with RF and vacuum

- 1.8.0mm& 13mm cartridge for fat removal;
- 2.body shaping vacuum rf body detoxing and shaping .



#### M5

- 1.Fractional and bipolar RF;
- 2.Plus and CTN operation mode;
3. intelligent tips recognize RF energy covering on skin fully or not;
- 4.24K golden pins two pins in one line,sensitive conduction is better.



#### Water and oxygen beauty equipment

- 1.spray jet :spray oxygen and nutrient on skin to promote absorption;
- 2.regular and constant injecting jet extract 99% unmixed oxygen and then inject it into skin ;
3. physicalogical saline solution jet:whiten skin color,moisture.



#### Qo

- 1.5 in 1 multi-function slimming beauty equipment;
- 2.adopt numerical rhythm control vacuum and ultrasound and RF and 905 soft laser and mechanical roller technolgy combined to beauty body ;
3. physical therapy.



#### BIO Micro-current

- 1.handle for face:activate cells,promote blood circulation,anti-inflammatory,conglutinant;
- 2.handle for hair:promote hair growth,relieve brain fatigue ;
- 3.handle for body:smoothing skin,lifting buttock,reducing wrinkle on hands and neck, firming and lifting breast.



#### Remage (thermage)

- 1.fractional and bipolar RF and cooling;
- 2.64pins & 144pins dot matrix handpieces ;
- 3.tightening and lifting face skin,brightening face skin color,restore the reticular layer of skin,wrinkles removal.



#### Fly Action

- 1.Vacuum & ultrasound & RF beauty equipment;
- 2.Slimming body,wrinkles removal,tightening skin,improve skin color & edema & microcirculation & stretch marks.



#### Dream Action Plus(portable)

- 1.Vacuum cavitation beauty equipment;
- 2.bipolar RF and numerical rhythm control vacuum technology combined;
3. physical therapy,contour body,beauty face& arms,dredge and metabolism.



#### Qo

- 1.905nm soft laser with vacuum handlepiece is good for neck detoxing,lymphatic drainage and body shaping;
- 2.40k Cavitation vacuum rf handlepiece is good for fat removal&body detoxing,shaping;
3. Mechanical roller vacuum rf handlepiece for body detoxing&shaping.



#### Water and oxygen beauty equipment

- 1.spray jet :spray oxygen and nutrient on skin to promote absorption;
- 2.regular and constant injecting jet extract 99% unmixed oxygen and then inject it into skin.